

# The new age in bioprocessing: Single-use bioreactors for microbial applications

Nico M.G. Oosterhuis<sup>1</sup>, Anton Tromper<sup>1</sup>, Nick A. van Biezen<sup>2</sup>

Compared to traditional bioreactors, single-use bioreactors show to be more flexible, diminish contamination risk, simplify validation and do require less infrastructure. However, the currently available single-use bioreactors are less suitable for application in microbial fermentation. The innovative CELL-tainer<sup>®</sup> bioreactor is multi-purpose and suited for both high-density mammalian cell cultures as well as microbial and fungal processes.

## Results

In the CELL-tainer<sup>®</sup> a significant improvement of the  $k_a$  value compared to the more traditional Wave Bioreactor<sup>™</sup> is found. When *E. coli* cells are cultivated at different scales and in different reactor types, it is shown that the CELL-tainer<sup>®</sup> bioreactor is well suited for *E. coli* cultivation. Optical densities measured in a 10L CELL-tainer<sup>®</sup> batch culture are comparable with a 1L and 100L stirred tank bioreactor. For all reactor types, air was enriched with oxygen. When no oxygen is added, growth is more slowly, but an OD = 70 is reached within 48 hrs in a fed-batch culture. The growth can be controlled by applying a gamma-irradiated in-line glucose measurement which is available as well.

## Costs

Due to the fast turn-around time and the restricted infrastructure needed, with the CELL-tainer<sup>®</sup> single-use bioreactor cost of operation can be reduced with at least 40% compared to an autoclavable or in-situ sterilizable fermenter.

As the system is very flexible in volume (0,2 – 15L in one-and-the-same bag) also for mammalian cultures cost of operation can be reduced even compared to the Wave type of single-use systems. Square shaped bags (3D) offer a volume range of 2,5L-25L in one-and-the-same bag for mammalian cell culture application and up to 20L for microbial fermentation.




## Conclusion

The CELL-tainer<sup>®</sup> single-use bioreactor is suitable for application with mammalian cell cultures as well as for microbial cultures. Applying in-line glucose measurement makes the system well suitable for process development. Also for screening purposes, pre-culture application and production of small batches the CELL-tainer<sup>®</sup> is very suitable, more flexible and reduces the risk of contamination. As the system also leads to significant cost reduction we expect a wide-spread of application in bio-processing and therefore a new age is started.

CELL-tainer<sup>®</sup> is a registered trademark of CELLution Biotech BV. Patents are pending.  
Data *E. coli* also courtesy to Lonza Hopkinton (as published at the IBC-Boston, 2009, WengLong Lin)

<sup>1</sup> CELLution Biotech BV, Dr. A.F. Philipsweg 15a, 9403AC Assen, The Netherlands,

<sup>2</sup> HAN BioCentre, P.O. box 6990, 6503GL Nijmegen, The Netherlands



The CELL-tainer<sup>®</sup> single use bioreactor is exhibited at STAND no. 4 (CELLON)



CELLution Biotech BV  
Dr. A.F. Philipsweg 15A  
9403AC Assen • The Netherlands  
T : +31 592 301929  
F : +31 592 308077  
info@cellutionbiotech.com • www.cellutionbiotech.com